

HOW TO SET UP YOUR KILOBASER

For further manuals & information please visit:
<https://manual.kilobaser.com>



GET STARTED

To watch our tutorial videos please visit:
<https://www.kilobaser.com/support>



Included in the package are:

- 1 Kilobaser device
- 1 Display
- 1 HDMI cable
- 1 USB-C to USB-A cable
- 1 Power supply with power cord
- 1 WIFI USB dongle
- 1 Blue 6mm tubing for gas supply 3m

ON / OFF Button

Press the button to turn on the device.

Touchscreen

Use the included HDMI and USB-C cable to connect the display to the device's back. (Connect the USB-C cable to the lower port on the display.)

ADDITIONAL QUESTIONS?

We are happy to assist with the setup of your specific application.

please contact:

support@kilobaser.com or
your local Kilobaser representative

SUPPORT CONTACT

phone: +43 660 3939 772

e-mail: support@kilobaser.com

web: www.kilobaser.com

address: Reininghausstrasse 13a
8020 Graz, Austria

